



Product Datasheet

FBN1 antibody (orb20385)


Description	Goat polyclonal antibody to FBN1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Target	FBN1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: Preliminary experiments in Human Brain (Cerebellum), Heart, Kidney, Skeletal Muscle and Skin lysates gave no specific signal but low background (at antibody concentration up to 1 µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	DASNIEDQSETEAN
Clonality	Polyclonal
MW	312
Entrez	2200
Dilution Range	ELISA: 1:2000, IF/ICC: 10µg/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U251 cell...



Immunofluorescence analysis of A431 cell...



Immunofluorescence analysis of Jurkat ce...