


Product Datasheet

CXXC4 antibody (orb20388)

Description	Goat polyclonal to CXXC4
Species/Host	Goat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	CXXC4
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: In paraffin embedded Human Breast shows cytoplasm staining in myoepithelial cells. Recommended concentration, 3-5µg/ml.WB: Preliminary experiments in Human Brain (Cerebellum, Hippocampus) lysates gave no specific signal but low background (at antibody concentration up to 1 µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	GTSLERTPVPSAEA
Clonality	Polyclonal
MW	21
Entrez	80319
Dilution Range	ELISA: 1:32000, WB: 1 µg/ml, IHC-P: 3-5µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...