

Product Datasheet

MDR1 antibody (Biotin) (orb20389)

Description	MDR1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Parasite
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Entamoeba histolytica Multidrug resistance protein 1 protein (1-114AA)
Target	Streptavidin
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P16875
Dilution Range	WB:1:500-5000, IHC-P:1:20-200, IF:1:50-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.