



Product Datasheet

KCNN2 antibody (orb20398)


Description	Goat polyclonal antibody to KCNN2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	KCNN2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 70kDa and 23kDa bands observed in Human Liver lysates (calculated MW of 63.8kDa according to NP_067627.2 and of 26.3kDa according to NP_740721.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	ESYDKHVTYNAER
Clonality	Polyclonal
MW	63.8; 26.3
Entrez	3781
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Liver lys...



Immunohistochemical staining of Human Th...



Immunohistochemical staining of Human Ki...