


Product Datasheet

ACADM antibody (orb20413)

Description	ACADM antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Medium-chain specific acyl-CoA dehydrogenase, mitochondrial protein (26-421AA)
Target	Pectate lyase 1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P11310
Dilution Range	ICC:1:10-100, IF:1:50-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human co...



Immunohistochemical staining of human ki...