

Product Datasheet

TAF8 antibody (orb20414)

Description

Goat polyclonal antibody to TAF8

Species/Host

Goat

Reactivity

Bovine, Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA

Applications

Target

TAF8 / TAFII43

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:16000. WB: Preliminary experiments gave an approx 38kDa band in Human Brain (Amygdala, Cerebellum, Hippocampus) lysates after 1 µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 34.3kDa according to NP_612639.2. The 38kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

KSYCEHTARTQPT

Clonality

Polyclonal

MW

34.3

Entrez

129685

Dilution Range

ELISA: 1:16000, WB: 1 µg/ml

Expiration Date

12 months from date of receipt.