


Product Datasheet

A4GNT antibody (orb20417)

Description	Goat polyclonal antibody to A4GNT
Species/Host	Goat
Reactivity	Bovine, Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	A4GNT / alpha4GNT
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 40kDa band observed in lysates of cell line K562 (calculated MW of 39.5kDa according to NP_057245.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	PRTYRDLIKGPEGS
Clonality	Polyclonal
MW	39.5
Entrez	51146
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb20417
(1ug/ml)
staining of
K562 lysat...