


Product Datasheet

CSMD1 antibody (orb20423)

Description	Goat polyclonal antibody to CSMD1
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CSMD1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000. WB: Approx 150kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 389kDa according to NP_150094.4). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Kraus et al, J Immunol. 2006 Apr 1;176(7):4419-30. PMID: 16547280). Recommended concentration: 1-3 µg/ml.
Protein Sequence	KGVREVNETHVTKTP
Clonality	Polyclonal
MW	389
Entrez	64478
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Cerebellu...