



Product Datasheet

SNX1 antibody (orb20425)


Description	Goat polyclonal antibody to SNX1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, WB
Target	SNX1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 75kDa band observed in Human lysates of cell line HeLa (calculated MW of 59.1kDa according to NP_003090.2). The observed molecular weight corresponds to earlier findings with different antibodies from other commercial sources. Recommended concentration: 1-3 µg/ml.
Protein Sequence	TNSSKPQPTYEELE
Clonality	Polyclonal
MW	59.1; 51.8; 53.3
Entrez	6642
Dilution Range	ELISA: 1:64000, WB: 0.1-1 µg/ml, IHC-P: 3-5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa lysate using SNX1 antibody.



Immunohistochemical staining of Human Uterus tissue using SNX1 antibody.



Immunofluorescence analysis of cells using SNX1 antibody.