


Product Datasheet

EYA1 antibody (orb20429)

Description	Goat polyclonal antibody to EYA1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	EYA1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 55kDa band observed in lysates of cell line HEK293 (calculated MW of 60.7kDa according to NP_742056.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	TDPTAEYSTIHSP
Clonality	Polyclonal
MW	61.2; 64.6; 60.7
Entrez	2138
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb20429
(2ug/ml)
staining of
HEK293 lys...