

Product Datasheet

SLC15A2 antibody (orb20435)

Description	Goat polyclonal antibody to SLC15A2
Species/Host	Goat
Reactivity	Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA
Target	PEPT2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Preliminary experiments in Human Brain (Frontal and Cerebral Cortex) and Kidney lysates gave no specific signal but low background (at antibody concentration up to 1 µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	QGNMIKLETKKTK
Clonality	Polyclonal
MW	81.8
Entrez	6565
Dilution Range	ELISA: 1:8000, WB: 1 µg/ml
Expiration Date	12 months from date of receipt.