



Product Datasheet

TGFBI antibody (orb20453)


Description	Goat polyclonal antibody to TGFBI
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Target	TGFBI
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000. WB: Approx 70kDa band observed in Human, Mouse and Rat Kidney lysates (calculated MW of 74.7kDa according to Human NP_000349.1). Recommended concentration: 0.2-0.5 µg/ml.
Protein Sequence	QLYTDRTEKLRPE
Clonality	Polyclonal
MW	74.7
Entrez	7045
Dilution Range	ELISA: 1:32000, WB: 0.2-0.5 µg/ml, IHC-P: 4-6 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Kidney (A...)



Western blot analysis of Human Kidney (A...)



Immunofluorescence analysis of paraformaldehyde-fixed cells (A...)