


Product Datasheet

LSP1 antibody (orb20457)

Description	Goat polyclonal to LSP1	
Species/Host	Goat	
Reactivity	Bovine, Canine, Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Target	LSP1 (aa261-275)	orb20457 (1ug/ml) staining of Peripheral...
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	This antibody is expected to recognize both reported isoforms (NP_002330.1; NP_001013271.1). Reported variants represent identical protein: NP_001013273.1, NP_001013271.1, NP_001013272.1 ELISA: Peptide ELISA: antibody detection limit dilution 1:16000. WB: Approx 55kDa band observed in Human Lymph Node, Tonsil and in Human Peripheral Blood Mononucleocytes lysates (calculated MW of 37.2kDa according to NP_002330.1). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Liu et al, J Exp Med. 2005 Feb 7;201(3):409-18. PMID: 15684321). Recommended concentration: 1-3 µg/ml.	
Protein Sequence	ASTKSRWETGEVQAQ	
Clonality	Polyclonal	
MW	37.2; 30.2	
Entrez	4046	
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml	
Expiration Date	12 months from date of receipt.	