


Product Datasheet

STAB1 antibody (orb20458)

Description	Goat polyclonal antibody to STAB1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	Stabilin 1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as large as the calculated size of 275kDa according to NP_055951.2. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring
Protein Sequence	DFPDTQRILTVK
Clonality	Polyclonal
MW	275
Entrez	23166
Dilution Range	ELISA: 1:8000, IHC-P: 3-5µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...