


Product Datasheet

RSL1D1 antibody (orb20472)

Description	Goat polyclonal antibody to RSL1D1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	RSL1D1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 75kDa band observed in lysates of cell line K562 (calculated MW of 55.0kDa according to NP_056474.2). The observed molecular weight corresponds to earlier findings with different antibodies from other commercial sources. Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	PTARKQLDKEQVRK
Clonality	Polyclonal
MW	55
Entrez	26156
Dilution Range	ELISA: 1:16000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb20472
(1µg/ml)
staining of
K562 lysat...