

Product Datasheet

TRIM29 antibody (orb20473)

Description

Goat polyclonal antibody to TRIM29

Species/Host

Goat

Reactivity

Bovine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, FC, IF, IHC, WB

Applications

Target

TRIM29 / ATDC

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:128000. WB: Approx 70kDa band observed in Human Peripheral Blood Mononucleocytes lysates (calculated MW of 65.8kDa according to NP_036233.2). Recommended concentration: 0.3-1 µg/ml. An additional band of unknown identity was also consistently observed at 37kDa. This band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any such results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

EKDRIKSFTTNE

Clonality

Polyclonal

MW

65.8; 37.3

Entrez

23650

Dilution Range

ELISA: 1:128000, WB: 0.1-1 µg/ml, IHC-P: 3.75 µg/ml


Expiration Date

12 months from date of receipt.

250kDa
150kDa
100kDa

50kDa
37kDa
25kDa
20kDa
15kDa

Western blot analysis
of Peripheral Bloo...



Immunohistochemical
staining of Human
Sk...