

Product Datasheet

IL5 Recombinant Protein (Pig) (orb2050547)

Description	IL5 Recombinant Protein (Pig)
Preservatives	Liquid or Lyophilized powder
Form/Appearance	Liquid or Lyophilized powder
Storage	-20°C or -80°C
Tag	N-terminal 6xHis-tagged
Note	For research use only
Protein Sequence	VENTMNRLVAETLTLLSIHRTLLIGDGNLMISTPVHTNHQLCIEEVFQGIDTLKNQ TARGDAVEKLFQNLSLIKEYIDRQKKNCGGERWRVTQFLDYLVFLGVINTEWT MES
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	17 kDa
Uniprot ID	Q9MYM5
NCBI	NP_999370
Expiration Date	6 months from date of receipt.