


Product Datasheet

DCUN1D1 antibody (orb20516)

Description	Goat polyclonal antibody to DCUN1D1
Species/Host	Goat
Reactivity	Bovine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	DCUN1D1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000. WB: Approx 28kDa band observed in HepG2 lysates (calculated MW of 30.1kDa according to NP_065691.2). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	RPQIAGTKSTT
Clonality	Polyclonal
MW	30.1
Entrez	54165
Dilution Range	ELISA: 1:16000, WB: 0.3-1 µg/ml, IHC-P: 5-10µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 lysate us...



Immunohistochemical staining of Human Ki...