

Product Datasheet

CPLX2 Recombinant Protein (Human) (orb2051815)

Description	CPLX2 Recombinant Protein (Human)
Preservatives	Liquid or Lyophilized powder
Form/Appearance	Liquid or Lyophilized powder
Storage	-20°C or -80°C
Tag	N-terminal GST-tagged
Note	For research use only
Protein Sequence	MDFVMKQALGGATKDMGKMLGGEEKDPDAQKKEEERQEALRQQEEERKAKH ARMEAEREKVRQQIRDKYGLKKKEEKEAEEKAALQPCESLTPKKAIPAGCGD EEEEEEESILDTVLKYLPGLQDMFKK
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	42.4 kDa
Uniprot ID	Q6PUV4
NCBI	NP_001008221
Expiration Date	6 months from date of receipt.