

Product Datasheet

ERAP1 Recombinant Protein (Mouse) (orb2052141)

Description	ERAP1 Recombinant Protein (Mouse)
Preservatives	Liquid or Lyophilized powder
Form/Appearance	Liquid or Lyophilized powder
Storage	-20°C or -80°C
Tag	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Note	For research use only
Protein Sequence	PCVQRAERYFREWKSSNGNMSIPIDVTLAVFAVGAQNTGWDFLYSKYQSSLSS TEKSQIEFSLCTSKDPEKLQWLLDQSFKEIIKTQEFPHILTLIGRNPVGYPLAWKF LRENWNKLQKFELGSSSIAHVMGTDDQFSTRARLEEVKGFFSSLKENGSQLR CVQQTETIEENIRWMDKNFDKIRLWLQKEKPELL
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	28.3 kDa
Uniprot ID	Q9EQH2
NCBI	NP_109636
Expiration Date	6 months from date of receipt.