

Product Datasheet

SULT1A1 Recombinant Protein (Rat) (orb2052180)

Description	SULT1A1 Recombinant Protein (Rat)
Preservatives	Liquid or Lyophilized powder
Form/Appearance	Liquid or Lyophilized powder
Storage	-20°C or -80°C
Tag	N-terminal 10XHis-SUMO-tagged and C-terminal Myc-tagged
Note	For research use only
Protein Sequence	MEFSRPPLVHVKGIPLIKFAETIGPLQNFTAWPDDLLISTYPKSGTTWMSEILDMI YQGGKLEKCGRAPIYARVPFLEFKCPGVPSGLETLEETAPRLLKTHLPLSLLPQSL LDQKVKVIYIARNAKDVVSYNMFYNMAKLHPDPGTWDSFLENFMDGEVSYGS WYQHVKWELRHHPVLYLFYEDIKENPKREIKKILEFLGRSLPEETVDSIVHHT SFKKMKENCMTNYTTIPTTEIMDHNVSFMRKGTGDWKNFTVAQNERFDAHY AKTMTDCDFKFRCEL
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	53.9 kDa
Uniprot ID	P17988
Expiration Date	6 months from date of receipt.