

Product Datasheet

ARMS2 antibody (orb20523)

Description

Goat polyclonal antibody to ARMS2

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC

Applications
Target

ARMS2

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:4000. WB: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as small as the calculated size of 11.4kDa according to NP_001093137.1. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

EGASDKQRSK

Clonality

Polyclonal

MW

11.4

Entrez
387715
Dilution Range


ELISA: 1:4000, IHC-P: 2-4 µg/ml

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of Human Re...



Immunohistochemical staining of Human Re...