

Product Datasheet

THY1 Recombinant Protein (Mouse) (orb2053190)

Description	THY1 Recombinant Protein (Mouse)
Preservatives	Liquid or Lyophilized powder
Form/Appearance	Liquid or Lyophilized powder
Storage	-20°C or -80°C
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Protein Sequence	QKVTSLTACLNVNQLRLDCRHENNTKDNSIQHEFSLTREKRKHVLSGTLGIPEHT YRSRVTLNQPVIKVLTLANFTTKDEGDYFCELQVSGANPMSSNKSISVYRDKLV KC
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	28.8 kDa
Uniprot ID	P01831
NCBI	NP_033408
Expiration Date	6 months from date of receipt.