



Product Datasheet

ANTXR2 antibody (orb20537)

Description	Goat polyclonal antibody to ANTXR2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	ANTXR2 / CMG2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 60kDa band observed in Human, Mouse and Rat Kidney lysates (calculated MW of approx 54kDa according to NP_477520.2; NP_001139266.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	HEGLKLANEQIKK
Clonality	Polyclonal
MW	53.5; 53.7
Entrez	118429
Dilution Range	ELISA: 1:16000, IHC-P: 2.5 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Liver tissue



Immunohistochemical staining of Human Kidney tissue