


Product Datasheet

MLLT10 antibody (orb20538)

Description	Goat polyclonal antibody to MLLT10
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	MLLT10 / AF10
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Experiments gave bands at approx 140kDa and 60kDa in Rat Testis lysates. These bands correspond to earlier findings with different antibodies from other commercial sources.. This protein has a calculated MW of 109kDa according to Human NP_004632.1, 107kDa according to Human NP_001009569.1and 103kDa according to Rat NP_001012162.1. Recommended concentration: 1-3 µg/ml.
Protein Sequence	RLEDTTARFTNAN
Clonality	Polyclonal
MW	109; 107
Entrez	8028
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb20538
(1ug/ml)
staining of
Rat Testis...