


Product Datasheet

GABRG2 antibody (orb20540)

Description	Goat polyclonal antibody to GABRG2
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	GABRG2 (aa400-410)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 50kDa band observed in Mouse and Rat Brain lysates (calculated MW of 55.1kDa according to Mouse NP_032099.1 and 54.1kDa according to Rat NP_899156.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	NNATHLQERDE
Clonality	Polyclonal
MW	55.2; 54.2; 59.3
Entrez	2566
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb20540
(2ug/ml)
staining of
Mouse Brai...