


Product Datasheet

TTC8 antibody (orb20541)

Description	Goat polyclonal antibody to TTC8
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	TTC8
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 60kDa band observed in Human Testis lysates (calculated MW of 58.4kDa according to NP_653197.2;). Recommended concentration: 0.2-0.6 µg/ml.
Protein Sequence	KEVLKQDNTHVE
Clonality	Polyclonal
MW	58.4; 57.3; 54.0
Entrez	123016
Dilution Range	ELISA: 1:8000, WB: 0.2-0.6 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Testis ly...