

Product Datasheet

NIPBL antibody (orb20548)

Description

Goat polyclonal antibody to NIPBL

Species/Host

Goat

Reactivity

Bovine, Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IF, IHC

Applications
Target

NIPBL (aa889-902)

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:32000. WB: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as large as the calculated size of 316kDa according to NP_597677.2 or 304kDa according to NP_056199.2. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

RPDSPRVKQGDSNK

Clonality

Polyclonal

MW

316; 304

Entrez


25836

Dilution Range


ELISA: 1:32000, IHC-P: 2 µg/ml

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of Human Ki...



Immunohistochemical staining of Human Ce...