

Product Datasheet

ACPP antibody (orb20559)

Description

Goat polyclonal to ACP

Species/Host

Goat

Reactivity

Bovine, Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IF, IHC, WB

Applications

Target

ACPP / PAP

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

IHC: In paraffin embedded Human Prostate shows cytoplasm staining in some secretory cells of the gland. Recommended concentration, 3-5µg/ml.WB: Approx 50kDa band observed in Human Prostate lysates (calculated MW of 44.6kDa according to NP_001090.2). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Patel et al, ,Clin Chem. 1986 Oct;32(10):1832-5. PMID: 3530530). Recommended concentration: 0.3-1 µg/ml.

Protein Sequence

YGIHKQKEKSRLQ

Clonality

Polyclonal

MW

44.6

Entrez

55

Dilution Range

ELISA: 1:64000, WB: 0.3-1 µg/ml, IHC-P: 3-5µg/ml


Expiration Date

12 months from date of receipt.

250kDa
150kDa
100kDa

50kDa
37kDa
25kDa
20kDa
15kDa

Western blot analysis
of human prostate ...



Immunohistochemical
analysis of formalin...