



Product Datasheet

ALDH3A1 antibody (orb20567)

Description	Goat polyclonal to ALDH3A1
Species/Host	Goat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	ALDH3A1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Reported variants represent identical protein: NP_001128639.1, NP_000682.3, NP_001128640.1 ELISA: Peptide ELISA: antibody detection limit dilution 1:32000. WB: Approx 55kDa band observed in Mouse Eye lysates (calculated MW of 50.5kDa according to Mouse NP_031462.2). Recommended concentration: 0.01-0.03 µg/ml.
Protein Sequence	RYPPSPAKMTQH
Clonality	Polyclonal
MW	50.4
Entrez	218
Dilution Range	ELISA: 1:32000, WB: 0.01-0.03 µg/ml, IHC-P: 3-5µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Eye lysat...



Immunohistochemical staining of Human Lu...