



## Product Datasheet

### MEIS2 antibody (orb20570)

|                            |                                                                                                                                                                                                                                                                                                                                                        |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Goat polyclonal antibody to MEIS2                                                                                                                                                                                                                                                                                                                      |
| <b>Species/Host</b>        | Goat                                                                                                                                                                                                                                                                                                                                                   |
| <b>Reactivity</b>          | Bovine, Canine, Human, Mouse                                                                                                                                                                                                                                                                                                                           |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                                                                                                                                                           |
| <b>Tested Applications</b> | ELISA, WB                                                                                                                                                                                                                                                                                                                                              |
| <b>Target</b>              | MEIS2                                                                                                                                                                                                                                                                                                                                                  |
| <b>Preservatives</b>       | Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.                                                                                                                                                                                            |
| <b>Storage</b>             | Aliquot and store at -20°C. Minimize freezing and thawing.                                                                                                                                                                                                                                                                                             |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                                                                                                                  |
| <b>Application notes</b>   | ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 50 + 45kDa bands observed in Rat Thymus lysates and in lysates of cell line KELLY (calculated MW of 43.8kDa according to Human NP_733777.1, 51.8kDa according to Human NP_733775.1 and of 51.0kDa according to Rat NP_001101228.1). Recommended concentration: 0.3-1 µg/ml. |
| <b>Protein Sequence</b>    | QSNRAGFLDPSVSQ                                                                                                                                                                                                                                                                                                                                         |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                                                                                                                                             |
| <b>MW</b>                  | 43.8;; 51.8; 42.3                                                                                                                                                                                                                                                                                                                                      |
| <b>Entrez</b>              | <a href="#">4212</a>                                                                                                                                                                                                                                                                                                                                   |
| <b>Dilution Range</b>      | ELISA: 1:16000, WB: 0.1-0.3 µg/ml                                                                                                                                                                                                                                                                                                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                                                                                                                                        |



Western blot analysis of Rat Thymus lysa...



Western blot analysis of Mouse Brain (La...