


Product Datasheet

Myogenin antibody (orb20580)

Description	Goat polyclonal to Myogenin
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Myogenin / MYF4
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 38kDa band observed in Mouse Skeletal Muscle lysates (calculated MW of 25.0kDa according to Human NP_002470.2 and of 25.2kDa according to Mouse NP_112466.1). The observed molecular weight corresponds to earlier findings with different antibodies from other commercial sources. Recommended concentration: 2-4 µg/ml.
Protein Sequence	TDAHNLHSLTSIVD
Clonality	Polyclonal
MW	25
Entrez	4656
Dilution Range	ELISA: 1:16000, WB: 2-4 µg/ml
Expiration Date	12 months from date of receipt.



orb20580
(2ug/ml)
staining of
Mouse Skel...