



Product Datasheet

ALDH3A2 antibody (orb20585)

Description	Goat polyclonal to ALDH3A2
Species/Host	Goat
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	ALDH3A2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize both reported isoforms (NP_001026976.1; NP_000373.1).ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 55kDa band observed in Mouse and Rat Liver lysates (calculated MW of 54.0kDa according to Mouse NP_031463.2 and 54.1kDa according to Rat NP_113919.2). Recommended concentration: 0.03-0.1 µg/ml.
Protein Sequence	KGLKGESVNLRY
Clonality	Polyclonal
MW	54.1
Dilution Range	ELISA: 1:32000, WB: 0.03-0.1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat Liver lysat...



Immunohistochemical staining of Human Ad...