


Product Datasheet

GRIN3B antibody (orb20586)

Description	Goat polyclonal antibody to GRIN3B
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	GRIN3B
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 110kDa band observed in Human Brain (Cerebellum and Cerebral Cortex) lysates (calculated MW of 113kDa according to NP_619635.1). Recommended concentration: 0.5-1.5 µg/ml.
Protein Sequence	EAPPHSGRPGSQ
Clonality	Polyclonal
MW	113
Entrez	116444
Dilution Range	ELISA: 1:8000, WB: 0.5-1.5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Cerebellu...