


Product Datasheet

ALDH1B1 antibody (orb20589)

Description	Goat polyclonal antibody to ALDH1B1
Species/Host	Goat
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ALDH1B1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.WB: Approx 60kDa band observed in Human Liver lysates (calculated MW of 57.2kDa according to NP_000683.3). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	DKEQFERVLGYIQ
Clonality	Polyclonal
MW	57.2
Entrez	219
Dilution Range	ELISA: 1:4000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Liver lys...