

Product Datasheet

NLRX1 antibody (orb20605)

Description	Goat polyclonal to NLRX1
Species/Host	Goat
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	NLRX1 / NOD9
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize both reported isoforms (NP_078894.2; NP_733840.1).ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 100kDa band observed in Human Breast cancer lysates (calculated MW of 108kDa according to NP_078894.2 and 102kDa according to NP_733840.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	NFSGETLDSTDPSN
Clonality	Polyclonal
MW	108; 102
Entrez	79671
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb20605
(2ug/ml)
staining of
Human
Brea...