


Product Datasheet

IMP3 antibody (orb20607)

Description	Goat polyclonal to IMP3
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	IMP3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.WB: Approx 22kDa band observed in Human Liver lysates (calculated MW of 21.9kDa according to NP_060755.1). Recommended concentration: 0.05-0.2 µg/ml.
Protein Sequence	QRREDYTRYNQLSR
Clonality	Polyclonal
MW	21.9
Entrez	55272
Dilution Range	ELISA: 1:4000, WB: 0.05-0.2 µg/ml
Expiration Date	12 months from date of receipt.



orb20607
(0.05ug/ml)
staining of
Human L...