


Product Datasheet

SNAP23 antibody (orb20617)

Description	Goat polyclonal to SNAP23
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	SNAP23
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize both reported isoforms (NP_003816.2; NP_570710.1).ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 26kDa band observed in Human Placenta lysates (calculated MW of 23.4kDa according to NP_003816.2). The observed molecular weight corresponds to earlier findings with different antibodies from other commercial sources. Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	QIKRITDKAD
Clonality	Polyclonal
MW	23.4; 17.8
Entrez	8773
Dilution Range	ELISA: 1:8000, WB: 0.1-0.3 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



orb20617
(0.1ug/ml)
staining of
Human Pl...