


Product Datasheet

EFNB2 antibody (orb20621)

Description	Goat polyclonal to EFNB2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	EFNB2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 35kDa band observed in Human, Mouse and Rat Lung lysates (calculated MW of 36.9kDa according to NP_004084.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	CPKVDSKTVGQYE
Clonality	Polyclonal
MW	36.9
Entrez	1948
Dilution Range	ELISA: 1:16000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat Lung lysate...