


Product Datasheet

ARHGDIG antibody (orb20624)

Description	Goat polyclonal to ARHGDIG
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	ARHGDIG
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 25kDa band observed in lysates of cell line HeLa (calculated MW of 25.1kDa according to NP_001167.2). Recommended concentration: 0.5-2 µg/ml.
Protein Sequence	EAVPEYRAPGRK
Clonality	Polyclonal
MW	25.1
Entrez	398
Dilution Range	ELISA: 1:16000, WB: 0.5-2 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa lysate using ARHGDIG antibody (orb20624).