


Product Datasheet

FANCL antibody (orb20626)

Description	Goat polyclonal to FANCL
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	FANCL
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize both isoforms (NP_001108108.1; NP_060532.2).ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 45kDa band observed in Human Bone Marrow lysates (calculated MW of 43.4kDa according to NP_001108108.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	QNLKDVLEIDFP
Clonality	Polyclonal
MW	43.4; 42.9
Entrez	55120
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb20626
(1µg/ml)
staining of
Human
Bone...