

## Product Datasheet

### CD5L Antibody : Biotin (orb2067535)

**Description**

CD5L Antibody : Biotin

**Species/Host**

Rabbit

**Conjugation**

Biotin

**Tested Applications**

ELISA, IHC, WB

**Immunogen**

The immunogen for anti-AIM antibody: antibody was raised against a 16 amino acid peptide near the amino terminus of human AIM.

**Preservatives**

Liquid. Purified antibody supplied in 1x PBS buffer.

**Form/Appearance**

Liquid. Purified antibody supplied in 1x PBS buffer.

**Storage**

All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.

**Note**

For research use only

**Application notes**

Application Info: AIM antibody can be used for the detection of AIM by Western blot at 1 - 2 ug/ml. (Optimal dilution should be determined by user)

**Clonality**

Polyclonal

**Uniprot ID**
[O43866](#)
**NCBI**
[NP\\_005885.1](#), [XP\\_005245659.1](#)
**Expiration Date**

12 months from date of receipt.