

Product Datasheet

TLR5 Antibody : FITC (orb2067609)

Description

TLR5 Antibody : FITC

Species/Host

Rabbit

Conjugation

FITC

Tested Applications

ELISA, IHC, WB

Immunogen

The immunogen for anti-TLR5 antibody: antibody was raised against a peptide corresponding to 16 amino acids near the center of human TLR5.

Preservatives

Liquid. Purified antibody supplied in 1x PBS buffer.

Form/Appearance

Liquid. Purified antibody supplied in 1x PBS buffer.

Storage

All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.

Note

For research use only

Application notes

Application Info: TLR5 antibody can be used for detection of TLR5 by Western blot at 1 to 2 ug/ml. (Optimal dilution should be determined by user)

Clonality

Polyclonal

Uniprot ID

[O60602](#)

NCBI

[NP_003259](#)

Expiration Date

12 months from date of receipt.