

Product Datasheet

PTK2 Antibody (Phospho-Tyr397) : HRP (orb2069225)

Description

PTK2 Antibody (Phospho-Tyr397) : HRP

Species/Host

Rabbit

Conjugation

HRP

Tested

IF, IHC-P, WB

Applications
Immunogen

KLH conjugated synthetic phosphopeptide derived from human FAK 1 around the phosphorylation site of Tyr397

Preservatives

Liquid. Purified antibody is supplied in high phosphate PBS, 100 mM phosphate, 150 mM NaCl, pH 7.6.

Form/Appearance

Liquid. Purified antibody is supplied in high phosphate PBS, 100 mM phosphate, 150 mM NaCl, pH 7.6.

Storage

All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.

Note

For research use only

Application notes

Application Info: WB(1:100-1000)IHC-P(1:100-500)IF(1:50-200)

Isotype

IgG

Clonality

Polyclonal

Uniprot ID
Q05397
Expiration Date

12 months from date of receipt.