

Product Datasheet

CXCL9 Antibody : HRP (orb2069444)

Description

CXCL9 Antibody : HRP

Species/Host

Rabbit

Conjugation

HRP

Tested

IF, IHC

Applications
Immunogen

This CXCL9 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 30-80 amino acids derived from human CXCL9.

Preservatives

Liquid. Purified antibody is supplied in high phosphate PBS, 100 mM phosphate, 150 mM NaCl, pH 7.6.

Form/Appearance

Liquid. Purified antibody is supplied in high phosphate PBS, 100 mM phosphate, 150 mM NaCl, pH 7.6.

Storage

All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.

Note

For research use only

Application notes

Application Info: IHC-P(1:100-500)IF(1:50-200)

Isotype

IgG

Clonality

Polyclonal

Uniprot ID
Q07325
Expiration Date

12 months from date of receipt.