

Product Datasheet

MC5R Antibody : HRP (orb2070044)

Description	MC5R Antibody : HRP
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	FC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human MC5 Receptor
Preservatives	Liquid. Purified antibody is supplied in high phosphate PBS, 100 mM phosphate, 150 mM NaCl, pH 7.6.
Form/Appearance	Liquid. Purified antibody is supplied in high phosphate PBS, 100 mM phosphate, 150 mM NaCl, pH 7.6.
Storage	All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4°C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.
Note	For research use only
Application notes	Application Info: WB(1:100-1000)FC(1:20-100)IHC-P(1:100-500)IF(1:50-200)
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.