

Product Datasheet

IRS2 Antibody : HRP (orb2070443)

Description

IRS2 Antibody : HRP

Species/Host

Rabbit

Conjugation

HRP

Tested

IF, IHC-P, WB

Applications
Immunogen

This IRS2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 540-590 amino acids from of human IRS2.

Preservatives

Liquid. Purified antibody is supplied in high phosphate PBS, 100 mM phosphate, 150 mM NaCl, pH 7.6.

Form/Appearance

Liquid. Purified antibody is supplied in high phosphate PBS, 100 mM phosphate, 150 mM NaCl, pH 7.6.

Storage

All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4°C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.

Note

For research use only

Application notes

Application Info: WB(1:100-1000)IHC-P(1:100-500)IF(1:50-200)

Isotype

IgG

Clonality

Polyclonal

Uniprot ID
[Q9Y4H2](#)
Expiration Date

12 months from date of receipt.