

Product Datasheet

CTDSP1 Antibody : HRP (orb2073113)

Description

CTDSP1 Antibody : HRP

Species/Host

Rabbit

Conjugation

HRP

Tested

ELISA, IHC, WB

Applications
Immunogen

The antiserum was produced against synthesized peptide between 22-71 amino acids from the N-Terminal region of human CTDSP1.

Preservatives

Liquid. Purified antibody is supplied in high phosphate PBS, 100 mM phosphate, 150 mM NaCl, pH 7.6.

Form/Appearance

Liquid. Purified antibody is supplied in high phosphate PBS, 100 mM phosphate, 150 mM NaCl, pH 7.6.

Storage

All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4°C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.

Note

For research use only

Application notes

Application Info: WB: 1:500 - 1:1000 IHC: 1:50 - 1:100 ELISA: 1:10000

Clonality

Polyclonal

MW

29 kDa

Uniprot ID
[Q9GZU7](https://www.uniprot.org/entry/Q9GZU7)
Expiration Date

12 months from date of receipt.