

Product Datasheet

BAZ2A Antibody : FITC (orb2073246)

Description	BAZ2A Antibody : FITC
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide between 1244-1293 amino acids from the Internal region of human BAZ2A.
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer.
Storage	All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.
Note	For research use only
Application notes	Application Info: WB: 1:500 - 1:1000ELISA: 1:10000
Clonality	Polyclonal
MW	208 kDa
Uniprot ID	Q9UIF9
Expiration Date	12 months from date of receipt.