

Product Datasheet

Nucleoside Phosphorylase antibody (orb20735)

Description

Rabbit polyclonal antibody to Nucleoside

Species/Host

Rabbit

Conjugation

Unconjugated

Preservatives

Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. no preservative added, as it may interfere with the antibody activity. No foreign protein added.

Form/Appearance

Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. no preservative added, as it may interfere with the antibody activity. No foreign protein added.

Storage

Storage: The lyophilised PAb fraction is shipped at ambient temperature and may be stored at +4°C; prolonged storage at or below -20°C. Prior to use, an aliquot is thawed slowly at ambient temperature, spun down again and used to prepare working dilutions by adding sterile phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided. Working dilutions should be stored at +4°C, not refrozen, and preferably used the same day. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance of the product

Note

For research use only

Clonality

Polyclonal

Source

Nucleoside phosphorylase isolated and purified from calf spleen. Freund's complete adjuvant is used in the first step of the immunization procedure.

Uniprot ID
P55859
Hazard
Information

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.